# EN User Manual

## LED face mask light therapy.

Model: MML

Thank you for purchasing our product. Please read the user manual carefully before first use. You should keep the manual for later reference as it contains important information.

### Package contents

1. LED mask
2. USB-A - USB-C power cable
3. Mounting band
4. User Manual

Check if the product is not damaged. Damaged equipment can be taken back to the point of purchase for repair or replacement in the case of damage and defects attributable to the product itself. For environmental protection purposes, dispose of the packaging in appropriate recycling containers. The packaging should not come into the hands of children. For the safety of children, please do not leave freely any accessible parts of the packaging (plastic bags, cardboard boxes, polystyrene, etc.). Risk of suffocation!

### Intended use

The device is used for facial light therapy and should only be used for cosmetic purposes, i.e. skin care. The product is intended for domestic, indoor use and cannot be used for professional purposes. The product is not waterproof and cannot be used while bathing or showering.

### Safety precautions

- Read the instructions carefully before use. Non-compliance with these instructions may be dangerous or illegal.
- The device should only be used for its intended purpose.
- The product is designed for facial skin care. It is not intended for laboratory or industrial use. It is not a medical device.
- Application time on the face must not exceed 20 minutes (for one session)
- Do not use in people with heart disease, patients with pacemakers or other electronic medical devices, people who drink alcohol and children
- The packaging is not a toy. After unpacking the product, keep the packaging away from the reach of children.
- The device is not a toy, make sure that it is not within the reach of children and they cannot play with it. The product is intended for adults.
- This product may only be repaired by a qualified manufacturer's service or an authorised service centre. Repair of the device by an unqualified service or by an unauthorised service centre may result in damage to the device and loss of warranty.
- The manufacturer is not responsible for any damage caused by self-modifications.
- Never use the product if it is damaged.
- Protect the device from damage due to falling and strong shocks.
- Prevent water from getting into the device. Do not immerse in water or other liquids. Do not use near sinks, bathtubs or other containers filled with water.
- Do not use on damaged, irritated skin, with active inflammatory lesions, herpes, after fresh treatments (acids, laser, micro-needling)
- If an allergic skin reaction occurs during or after using the product, discontinue use and consult a doctor.
- Use the device at room temperature.
- Wash and dry hands before use.
- Do not use on wet hair
- Do not look directly at the emitted light
- Don't sleep in a mask, don't use while charging
- The effect of the device may vary depending on skin type, extent and frequency of use.
- Do not use in case of burns, post-traumatic infections or other skin problems.
- The product is not suitable for pregnant women.
- For hygienic reasons, the device should only be used by one person.
- When charging the battery, connect the device in such a way that the cable cannot be pinched, doesn't hang, or lie on sharp edges.
- Plug the device into an easily accessible socket so that it can be quickly disconnected from the power supply in the event of an emergency.
- The AC adapter is not included, connect the power cable to the sources according to the parameters specified in the specification.
- In the event of disturbing noise or smoke coming from the device, unplug it immediately.
- Use a dry and soft cloth or a slightly damp cloth (not soaked in chlorine or solvents) to clean the device.
- Do not use highly acidic or alkaline agents.

### Battery safety precautions:

- The device is powered by a built-in rechargeable battery, integral to the device, and is not within the scope of replacement by the user.
- Avoid exposing batteries to extremely low or high temperatures (below 10°C / 50°F or above 45°C / 113°F). Extreme temperatures can affect the battery life span.
- Keep batteries away from liquids and metal objects as this may lead to complete or partial damage.
- Use the battery only for its intended purpose. Do not destroy, damage or throw batteries into the fire – this can be dangerous and may cause fire. Explosion hazard!
- Ingesting batteries may be fatal! Keep them away from children and pets. Seek medical attention immediately if you swallow a battery.
- Do not dispose of used or damaged batteries with household waste. They must be handed over to a designated recycling point in the local municipality.

- Do not dismantle, modify or short-circuit the battery terminals!
- Before first use, charge the battery in accordance with the operating instructions.

**NOTE:** If you do not use the device for an extended period, charge it at least once every 3 months nonetheless. Do not allow the battery to discharge completely.

**NOTE:** The battery is non-replaceable and is an integral part of the device. Do not attempt to replace it.

## Construction

### 1. FUNCTION BUTTON

### 2. POWER SOCKET

A function button allows you to switch the mask on and off and also change the colour of the light.

## Technical specification

battery:	Li-poly 3.7V --- 400 mAh 1.48 Wh
Power supply:	5V --- 1A
charging port :	USB-C charging port
Number of light modes:	7 light modes
Dimensions:	22.8 cm x 18.3 cm x 13.5 cm

## Types of lights and their operation

**RED LIGHT:** stimulates fibroblasts to produce collagen and elastin, promotes healing, reduces fine lines and improves blood circulation and skin radiance.

**GREEN LIGHT:** soothes the skin, helps even out skin tone and gently promotes the reduction of post-inflammatory hyperpigmentation. **BLUE LIGHT:** has an anti-acne effect by acting on *C. acnes* bacteria, reducing inflammation and the number of blemishes. **VIOLET/BLUE-VIOLET LIGHT:** a combination of red and blue, it combines an anti-acne effect with gentle regeneration and support for the brightening of acne marks. **LIGHT BLUE (CYAN):** soothes irritation and redness, promotes skin comfort after treatments. **LIGHT GREEN:** a variant of green with a soothing, toning and unifying effect. **WHITE LIGHT:** broad spectrum (all colours), general mode for brightening and balancing the appearance of the skin.

## Recommended combinations

For active acne - blue followed by (optionally) purple or a short session of red; for anti-ageing and for glow - red alone or red followed by white; for hyperpigmentation and uneven skin tone - green/light green alternating with red (e.g. every other day); for sensitive skin or after treatments - light blue (cyan) or green in shorter sessions.

## Safety instructions:

Use on clean, dry skin; session 10-20 minutes, 3-5 times a week for at least 6-8 weeks; protect eyes, especially under blue and violet light; apply moisturising/antioxidant serum when finished; discontinue if irritated; avoid while taking photosensitising medication and consult use with your doctor if you have a skin/eye condition.

## Charging the device battery

Before first use, charge the device's battery. Connect the USB-A plug of the power cable to a power source with the specifications indicated, 5V ---. Plug the other end of the cable with a USB-C plug into the power socket on the device. Then plug the charger into a mains socket, the charging indicator, the LED located under the function button will light up red (the LED lights up very softly). When fully charged, the indicator light under the function button will go out. Install the mounting band on all 4 brackets on the mask.

**NOTE! The charger is not included.**

## Use

To switch the device on or off, press and hold the function key. Once the device is switched on, a blue LED will light up under the button and the LED mask will start in red light mode. Changing the LED light is possible by pressing the function button. The light changes sequentially: red, green, blue, violet, light green, blue, white. The mask is recommended not to be used for more than 20 minutes at a time or more than once a day. Prepare your face, apply make-up remover, cleanse and dry your face, you can apply a light serum don't use with thick, oil-based cosmetics. Remove any jewellery from your face. Apply the mask, adjust the straps so that the mask fits firmly but without pressure. Close your eyes, switch on the device and select the appropriate mode.

**Don't sleep in a mask and don't use while charging.**

## Maintenance and storage

Clean the device thoroughly after each use. Turn off the appliance. Use a dry and soft cloth or a slightly damp cloth (not soaked in chlorine or solvents) to clean the device. The device is not waterproof and cannot be submerged or rinsed with water. Do not clean with caustic chemicals. Observe general health and safety regulations and exercise particular caution when cleaning. Do not perform maintenance while the device is charging. Store the device in a dry, ventilated place, away from heat sources. If you do not intend to use the products for a long time, it is advisable to store them in their original packaging.

## Disposal

The device is labelled with the symbol of a crossed-out rubbish bin, in accordance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE). The products marked with this symbol should not be disposed of with other household waste after their time of use. The user is obliged to dispose of waste electrical and electronic equipment by delivering it to the designated point where such hazardous

waste is recycled. Collecting this type of waste separately and ensuring it is properly processed contributes to the conservation of natural resources. Proper recycling of waste electrical and electronic equipment has a positive impact on human health and the environment. In order to obtain information on where and how to dispose of waste electronic equipment responsibly and in an environmentally safe manner, the user should contact the relevant local authority, the waste collection point or the point of sale where the equipment was purchased.

In accordance with EU Regulation 2023/1542 on batteries and waste batteries, this product is marked with a crossed-out rubbish bin symbol. The symbol indicates that the batteries used in this product should not be disposed of with normal household waste, but treated in accordance with the directive and local regulations. Batteries must not be thrown out with unsorted municipal waste. Users must use the available battery waste collection network that enables return, recycling, and disposal. In the EU, the collection and recycling of batteries is subject to separate procedures. To learn more about local battery recycling procedures, contact your municipality, waste management authority or landfill site.

## Explanation of symbols

The product complies with European technical standards defined as harmonised European Standards. The device complies with the directive 2014/30/EU. The declaration of conformity and the CE marking were issued.

The symbol means that the product must not be disposed of with other household waste. Any electrical equipment that is not recycled pollutes the air, soil and groundwater. The Polish law prohibits, under penalty of a fine, combining electrical and electronic equipment with other waste.

The symbol indicates the type of material used for the packaging - corrugated cardboard.

The symbol indicates the recyclability of the packaging.

The symbol indicates that you should read the instruction manual before use.

## Device manufacturer

### TelForceOne S.A.,

ul. Krakowska 119 Street, 50-428 Wrocław, Poland

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